







8165 Genu Neurexa

DE Gebrauchsanweisung	4
EN Instructions for use	8
FR Instructions d'utilisation	11
IT Istruzioni per l'uso	15
ES Instrucciones de uso	19
PT Manual de utilização	23
NL Gebruiksaanwijzing	27
SV Bruksanvisning	31
DA Brugsanvisning	34
NO Bruksanvisning	38
FI Käyttöohje	42
PL Instrukcja użytkowania	45
HU Használati utasítás	49
CS Návod k použití	53
SK Návod na používanie	57
EL Οδηγίες χρήσης	60
RU Руководство по применению	64
JA 取扱説明書	69

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INFORMATION

Date of last update: 2021-08-04

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident in connection with the product, in particular a worsening of the state of health, to the manufacturer and to the relevant authority in your country.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 8165 Genu Neurexa knee orthosis.

2 Intended use

2.1 Indications for use

The orthosis is intended **exclusively** for orthotic fittings of the lower limbs and **exclusively** for contact with intact skin.

The orthosis must be used in accordance with the indications.

2.2 Indications

- Paresis of leg muscles resulting in hyperextension of the knee

Caused for example by:

- Stroke
- Intervertebral disc prolapse in the lumbar spine
- Peripheral nerve damage

2.3 Contraindications

2.3.1 Absolute Contraindications

Not known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent scars that are swollen, reddening and hyperthermia of the fitted limb/body area; lymphatic flow disorders, including unclear soft tissue swelling distant to the body area to which the medical device will be applied; sensory and circulatory disorders in the leg area.

2.4 Effects

- Supports and stabilises the knee joint
- Improves proprioception
- Supports sensorimotor function
- Prevents hyperextension
- Promotes the reduction of oedemas and haematomas
- Can relieve pain

3 Safety

3.1 Explanation of warning symbols

Warning regarding possible risks of accident or injury.

Warning regarding possible technical damage.

3.2 General safety instructions

⚠ CAUTION

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- ▶ The product may be used by one person only.
- ▶ Clean the product regularly.

⚠ CAUTION

Mechanical damage to the product

Injuries due to changes in or loss of functionality

- ▶ Use caution when working with the product.
- ▶ Check the product for proper function and readiness for use.
- ▶ In case of changes in or loss of functionality, discontinue use of the product and have it checked by authorised, qualified personnel.

⚠ CAUTION

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage.

- ▶ Keep the product away from open flames, embers and other sources of heat.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilisation due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

NOTICE

Improper use and changes

Change in or loss of functionality as well as damage to the product

- ▶ Use this product with care and only for its intended purpose.
- ▶ Do not make any improper changes to the product.

4 Use

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel according to the instructions of the treating physician.
- ▶ Consult a physician if any exceptional changes are noted (such as worsening of the complaint).

4.1 Size selection

- 1) Measure the circumference of the lower leg **15 cm** below the centre of the patella (circumference A).
- 2) Measure the circumference of the thigh **15 cm** above the centre of the patella (circumference B).
- 3) Determine the size of the orthosis (see size chart).

4.2 Adaptation

NOTICE

Incorrect setting

Damage to the orthosis due to overloading of the material and improper fit of the orthosis due to breakage of load-bearing components

- ▶ The orthosis may only be adjusted by qualified personnel.
- ▶ Do not make any improper changes to the settings.

Adapting the splints

Individual adaptation of the aluminium bars is carried out by trained, qualified personnel.

- 1) Open the hook-and-loop cover for the bars (see fig. 1).
- 2) Take the bars out of the pockets (see fig. 2).
- 3) **INFORMATION: Do not bend the bars in the area of the joint; maintain a minimum distance of 3 mm from the plastic housing. Failure to comply could reduce the stability of the bars.** To optimise the fit, adapt the bars to the patient's leg contour using bending irons.
- 4) Reinsert the bars into the orthosis and fasten the hook-and-loop cover.

Adjusting the extension and/or flexion stops

- ▶ Change the extension and flexion stops only if advised to do so by the attending physician.

The markings on the stops indicate the angle at which the orthosis stops the leg in flexion or extension (see fig. 3).

- 1) Insert the required extension stop, alternatively the flexion stop (see fig. 4).
→ The stop will engage with slight resistance.
- 2) Move the orthotic joint and check the position of the stops.

4.3 Application

⚠ CAUTION

Incorrect or excessively tight application

Risk of local pressure and constriction of blood vessels and nerves due to improper or excessively tight application

- ▶ Ensure that the product is applied properly and fits correctly.

- 1) Apply the orthosis around the leg from behind and position the patella in the opening (see fig. 5).
- 2) Use the elasticity of the material and pull the protective edges over each other. Fasten the orthosis with the hook-and-loop closures on the lower leg, then repeat this process on the thigh (see fig. 5).
- 3) Pull the four fastening straps through the strap guide loops, tighten and close them (see fig. 6).
- 4) To adjust the extension limitation, loosen the crossed hyperextension strap in the knee fold (see fig. 7).
- 5) **INFORMATION: For maximum extension limitation, slightly flex the knee joint while fastening the crossed hook-and-loop closure.** Adjust the strap to the desired length and fasten the hook-and-loop closure (see fig. 7, see fig. 8).

4.4 Delivery

- 1) Check the product fits properly on the leg when handing the orthosis over to the patient.
- 2) Ask the patient to sit down, stand up and walk a few steps.
→ The centre of the orthotic joint should be at the height of the middle of the patella.
→ The orthosis should not be twisted on the leg.
→ The straps should be tight but not impair the patient's comfort.

Inform the patient:

Inspect the orthosis for damage

- ▶ Inspect the orthosis for damage daily. Contact the qualified personnel promptly if changes are noted.

5 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- ▶ Only clean the product with the approved cleaning agents.

Clean the orthosis regularly:

▶ **NOTICE! Detergent residues may cause skin irritation and wear of material.**

- 1) Take the bars out of the pockets.
- 2) Fasten all hook-and-loop closures.
- 3) Recommendation: Use a laundry bag or net.
- 4) Wash the orthosis in warm water at **30 °C** with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 5) Allow to air dry. Avoid exposure to direct heat (e.g. sunlight, oven or radiator).
- 6) Reinsert the bars into the orthosis.

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal information

All legal conditions are subject to the respective national laws of the country of use and may vary accordingly.

7.1 Liability

The manufacturer will only assume liability if the product is used in accordance with the descriptions and instructions provided in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorised modification of the product.

7.2 CE conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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