

# FORMFIT® TRACKER

The Formfit® Tracker with micro-adjustable ReflexWing™ buttress realigns the knee cap, providing dynamic lateral support for common knee cap conditions. Its FlyWeight™ fabric delivers lightweight breathability and a low-profile design that feels great against the skin and minimises migration. PowerLock™ straps work in tandem with heat mouldable CustomFit™ hinges to provide a customisable fit for your patient.

## INDICATIONS FOR USE

For knee conditions that may benefit from realignment of the knee cap and/or patients that would benefit from an added lateral support of the knee cap, such as:


- lateral knee cap subluxation or dislocation
- lateral maltracking

## FEATURES AND BENEFITS

- CustomFit™ hinges for customised fit and support
- Anti-Migration Tech™ minimises migration
- PowerLock™ straps for secure lateral support
- ReflexWing™ buttress for comfortable knee cap realignment
- FlyWeight™ fabric for soft feel and breathability
- Durable and machine washable



## REQUIRED MEASUREMENTS & ORDERING INFORMATION



Calf circumference 15 cm (6") below mid patella



Thigh circumference 15 cm (6") above mid patella

Size	Calf circumference 15 cm (6") below mid patella		Thigh circumference 15 cm (6") above mid patella		Left	Right
	cm	inch	cm	inch		
XS	25 - 30	10 - 12	33 - 39	13 - 15½	719162	719152
S	30 - 35	12 - 14	39 - 47	15½ - 18½	719163	719153
M	35 - 40	14 - 16	47 - 53	18½ - 21	719164	719154
L	40 - 46	16 - 18	53 - 60	21 - 23½	719165	719155
XL	46 - 51	18 - 20	60 - 67	23½ - 26½	719166	719156
XXL	51 - 56	20 - 22	67 - 75	26½ - 29½	719167	719157
3XL	56 - 61	22 - 24	75 - 81	29½ - 32	719168	719158

## ACCESSORIES

Description	Part#
ReflexWing™ Buttress	719100